




# blink<sup>®</sup>medical

the single-use instrument company



Text size is approximately  
2mm - raised up from the  
end face



DETAIL B  
SCALE 3 : 1

[www.blinkmedical.com](http://www.blinkmedical.com)



# Welcome

Blink Medical is one of the worlds largest manufacturers of quality, sterile, single use surgical instruments and bespoke instrument sets and packs. We supply instruments from our UK headquarters, to both NHS and private sector hospitals for use in either theatres or out patient departments. Through our exclusive distributor network, we export to over 20 countries worldwide including Australia, New Zealand and throughout continental Europe.

Offering the highest quality instruments, we provide products for Ophthalmic, Dermatology, Plastic Surgery, Cosmetic Surgery and ENT departments. More recently we have also introduced a wider range of Bi Polar Forceps and we have increased the range of standard off the shelf packs.

Since launching in 2004, our primary objective has been to supply the highest quality surgical instruments, in conjunction with first class customer service. Our high standards and attention to detail are reflected throughout the company, beginning with our manufacturing sites which adhere to ISO9001 and ISO13485, through to quality inspection, packaging, sterilisation, and the onward supply to the end user.

If you require any further information please call our customer service team on +44 121 386 8433. Alternatively to view our extensive product range visit us at [www.blinkmedical.com](http://www.blinkmedical.com) to see why we are known as...

*"the single use instrument company".*



**blink**medical  
the single-use instrument company





## Manufacturing Facility

All Blink surgical instruments are manufactured to the tightest tolerances. Each of our four main manufacturing sites are ISO compliant and visited regularly to ensure they work to Blinks demanding standards. Each instrument is inspected and ultrasonically cleaned before being shipped so that it arrives at our distribution centre in perfect condition.



## Logistics

Our UK Logistics Centre handles thousands of shipments every year, sending consignments to the UK and around the world. Whether ordering one box of forceps or a large consignment of customised sets, our highly trained professional staff are here to ensure your shipment is expertly handled and delivered on time.



## UK Quality Control

Our UK based quality control department is staffed by people with a wealth of experience in this vital area. Our two stage process includes multiple checks to ensure each instrument arrives in perfect condition. Gamma Radiation and Ethylene Oxide sterilisation are used to ensure complete effective sterility and cleanliness. Furthermore regular Bio Burden tests are also performed to ensure the correct level of sterilisation is achieved from one batch to another.



## Customer Service

Our customer services department will handle your order professionally and courteously, expediting it in the quickest time possible. Should you have a query or require a quotation, our friendly, helpful staff will aim to deal with your request politely and quickly. Where possible we will process any quotation for a bespoke set the same day.



## International Distribution

Our exclusive International Distributor network covers the majority of Europe, Australia and New Zealand. Our main distributors carry the Ophthalmic range of instruments, however some also now carry the Cosmetic, ENT and Dermatology lines as well. Should you represent a company interested in becoming a Blink Distributor or alternatively if you are a potential customer wishing to locate your nearest Blink distributor, please visit our website at [www.blinkmedical.com](http://www.blinkmedical.com) and click on the appropriate tab.



**10-1415 - Forcep Birks Suture Tying**  
8mm Tying Platform, O/L 110mm



**10-1302 - Forcep Utility Lister**  
1:2 Teeth 0.8mm with 6mm Tying Platform, O/L 114mm



**10-1203 - Forcep St Martins**  
1:2 Teeth 0.6mm, O/L 90mm  
**10-1204 - Forcep St Martins**  
1:2 Teeth 0.3mm, O/L 90mm



**10-1205 - Forcep Barraquer Fixation**  
9:10 Teeth, O/L 103mm



**10-1104 - Forcep Bonn**  
1:2 Teeth, 0.8mm with 6mm Tying Platform, O/L 82mm




**10-1105 - Forcep Bonn**  
1:2 Teeth, 0.12mm with 6mm Tying Platform, O/L 82mm

**10-9007 - Forcep Tying Mcpherson**  
Straight Plain Tip 6mm Tying Platform O/L 85mm



**10-9008 - Forcep Tying Mcpherson**  
Straight Plain Tip 6mm Tying Platform O/L 110mm



**10-8008 - Forcep Hoskins Straight Notched**  
0.5mm Notch 5mm Tying Platform O/L 90mm



**11-7009 - Forcep Pierse Collibri Notched**  
0.3mm Notch O/L 85mm



**11-7010 - Forcep Collibri Toothed**  
1:2 Teeth 0.8mm Tooth O/L 65mm  
**11-7011 - Forcep Collibri Toothed**  
1:2 Teeth 0.8mm Tooth O/L 110mm



**11-6010 - Forcep Tying Kelman McPherson**  
8mm Tying Platform O/L 85mm





### 11-5011 - Forcep Capsulorhexis Inamura

Flat Handle  
Cross Action  
45 Degree Shaft  
0.4mm Sharp Tip  
O/L 105mm



### 11-5012 - Forcep Capsulorhexis Utrata

Round Handle  
Direct Action  
45 Degree Shaft  
0.4mm Tip  
O/L 103mm



### 11-5013 - Forcep Capsulorhexis Utrata

Round Handle  
Direct Action  
Curved Shaft  
0.4mm Tip  
O/L 103mm



### 11-5020 - Forcep Capsulorhexis Inamura

Round Handle  
Cross Action  
Curved Shaft  
0.4mm Sharp Tip  
O/L 105mm



### 11-5021 - Forcep Capsulorhexis Inamura

Flat Handle  
Cross Action  
Curved Shaft  
0.4mm Sharp Tip  
OL 105mm



### 11-5022 - Forcep Capsulorhexis Utrata

Flat Handle  
Direct Action  
45 Degree Shaft  
0.4mm Sharp Tip  
OL 105mm



### 11-5023 - Forcep Capsulorhexis Micro Incision 1.8mm

Direct Action  
Curved Shaft  
0.4mm Sharp Tip  
OL 97mm




### 11-5024 - Forcep Capsulorhexis Utrata Curved

Flat Handle  
Direct Action  
Curved Shaft  
0.4mm Sharp Tip  
O/L 105mm



### 11-5025 - Forcep Capsulorhexis Fulcher

Flat Handle  
Direct Action  
Curved Shaft  
0.4mm Rounded Tip  
O/L 105mm






**10-1300 - Forcep Epilation**


Angled Tips  
O/L 91mm

New



**10-1301 - Forcep Suture Moorfields**

16mm Serrated Tying Platform  
O/L 107mm



**20-7001 - Scissor Iris Straight**

Rounded Tip O/L 100mm




**20-7002 - Scissor Iris Straight**

Sharp Tip O/L 100mm




**10-1311 - Forcep Suture Moorfields Utility**

14mm Serrated Tying Platform  
O/L 107mm



**11-1005 - Forcep Mosquito Curved**

Curved Tip  
O/L 100mm



**20-6003 - Scissor Vannas Straight**

O/L 95mm




**21-6004 - Scissor Vannas Curved**

O/L 95mm



**10-1006 - Forcep Mosquito Straight**

O/L 100mm



**20-4006 - Scissor Westcott Straight**

Rounded Tip O/L 120mm



**20-4007 -- Scissor Westcott Straight Sharp**

O/L 120mm





**21-4008 - Scissor Westcott Curved**  
Sharp Tip O/L 120mm



**21-4009 - Scissor Westcott Curved**  
Rounded Tip O/L 120mm




**21-4010 - Scissor Westcott Curved**  
Rounded Tip O/L 120mm



**21-5001 - Scissor Westcott Curved Small**  
Sharp Tip O/L 100mm




**21-3010 - Scissor Ong**  
Sharp Tip O/L 92mm




**40-5002 - Manipulator Lens Hook Sinskey Straight**  
Tip L 0.5mm O/L 124mm



**41-5001 - Manipulator Lens Hook Sinskey Angled**  
Tip Angled at 45 Degrees  
Tip L 0.5mm O/L 120mm



**41-3004 - Manipulator Spatula**  
Tip Angled at 45 Degrees  
Tip L 12mm W0.75mm O/L 135mm



**41-2005 - Rotator Nucleus Y Shaped**  
On horizontal axis  
Tip Angled at 45 Degrees  
Tip L 10mm O/L 124mm



**41-1006 - Lester IOL Manipulator**  
Tip Angled at 45 Degrees  
Tip L 12mm 1mm Branch O/L 125mm



**41-7008 - Chopper Greens Phaco**  
Tip Angled at 45 Degrees  
Tip L 10mm O/L 117mm




**41-7009 - Chopper Claw**  
Tip Angled at 45 Degrees  
Tip L 10mm O/L 123mm




**41-4005 - Chopper Sadiq**  
Tip Angled at 45 Degrees  
Tip L 11mm O/L 121mm



**41-4003 - Manipulator Phaco Chopper Mackool**  
Tip Angled at 45 Degrees  
Tip L 2mm O/L 120mm



**41-3005 - Neuhann Spatula**  
Tip Angled at 45 Degrees



**41-2006 - Rotator Nucleus Y Shaped**  
On vertical axis  
Tip Angled at 45 Degrees  
Tip L 10mm O/L 124mm



**41-6007 - Paddle Manipulator Drysdale**  
Tip Angled at 45 Degrees  
Tip L 11mm Spatula 1mm O/L 120mm




**41-8009 - Manipulator Mushroom Tipped**  
Tip Angled at 45 Degrees  
Tip L 11mm O/L 123mm



**41-4004 - Chopper Ambidextrous**





**41-4006 - Spatula Govan**  
Tip Angled at 45 Degrees  
Tip L 10mm O/L 130mm



**41-4007 - Haefliger Chopper Sadiq**  
Tip Angled at 45 Degrees  
Tip L 11mm O/L 115mm



**41-9010 - Manipulator Melonies Double Ended - Tip Angled at 45 Degrees**  
Sinskey Angled Tip L 0.5mm  
Mushroom Tip L 11mm O/L 123mm






**30-4000 - Speculum Barraquer**  
Solid Blade Angled Temporal L40mm W42mm




**30-4001 - Speculum Adjustable**  
Sliding Angled Temporal L45mm



**30-3001 - Speculum Burratto**  
Straight L57mm




**31-2003 - Speculum Barraquer**  
Angled Temporal L40mm W30mm



**30-2004 - Speculum Barraquer Straight**  
L40mm W30mm



**31-1006 - Speculum Kratz-Barraquer**  
Angled Temporal L40mm W30mm



**30-3002 - Speculum Paediatric Spring**  
Solid Blade 7mm O/L35mm L30mm W30mm




**31-3003 - Speculum Cook Adjustable**  
Solid Blades L40mm




**31-3004 - Speculum Clarkes**  
Angled Temporal L80mm




**30-1007 - Speculum Kratz-Barraquer**  
Straight L40mm W30mm




Temporal




Straight




**31-3100 - Speculum Neo Natal**  
Angled Temporal W21mm O/L34mm



**31-3101 - Speculum Neo Natal**  
Strong Spring Angled Temporal W21mm O/L34mm



**31-3002 - Speculum Leibinger**  
Angled Temporal O/L65mm






51-0000 - Flap Lifter Lasik - O/L122mm - Tip L12mm



51-0001 - Epithelial Spatula - O/L113mm - Spatula L7mm




51-0002 - Alcohol Well - O/L114mm I/D 9mm




The Blink range of Vitreoretinal Forceps have been designed to create a tactile feel to emulate that of a true reusable instrument. Competitively priced, the five instruments enable the surgeon to perform their procedure with confidence, knowing each instrument has been hand crafted, inspected and tested here in the UK.



60-1000 - Forcep End Gripping  
Tip 23g O/L132mm



60-1001 - Forcep ILM  
Tip 23g O/L132mm



60-1002 - Forcep A-Symetrical  
Tip 23g O/L132mm



60-1003 - Forcep Straight  
Tip 23g O/L132mm



60-2000 - Scissor Straight  
Tip 23g O/L132mm



**90-2008 - Probe Liebricht 000+00**  
O/L125mm



**90-2009 - Probe Liebricht 00+0**  
O/L125mm



**90-2010 - Probe Liebricht 1+2**  
O/L125mm



**90-2011 - Probe Liebricht 3+4**  
O/L125mm



**90-1011 - Dilator Punctal Double Ended**  
O/L138mm



**90-1012 - Nettleships Dilator**  
O/L116mm



**90-8012 - Bulldog Clamp**



**90-9013 - Desmares Retractor**  
14mm



**90-7001 - Needle Holder Castroviejo**  
Straight with lock O/L135mm 12mm Tying Platform



**90-7002 - Needle Holder Castroviejo**  
Without lock O/L135mm 12mm Tying Platform



**91-7003 - Needle Holder Castroviejo**  
Delicate Curve O/L135mm 12mm Tying Platform



**90-7004 - Needle Holder Troutman**  
O/L105mm 6mm Tying Platform



**90-5004 - Hook Stabismus Graefe Medium**  
O/L140mm Tip L10mm



**90-5005 - Hook Alexandra Muscle**  
O/L130mm Tip L11mm





**90-4005 - Clamp Chalazion Small**  
O/L90mm Clamp Diameter 11mm



**90-4006 - Clamp Chalazion Medium**  
O/L80mm Clamp Diameter 18mm x11mm



**90-4007 - Clamp Chalazion Large**  
O/L90mm Clamp Diameter 23mm x15mm



**90-1014 - Caliper and Rule**



**90-1015 - Double Ended Marker**  
3.5 - 4.0mm O/L 113mm



**90-4008 - Desmares Entropian Clamp Left**  
O/L95mm Clamp Diameter 24mm x 16mm



**90-4009 - Desmares Entropian Clamp Right**  
O/L95mm Clamp Diameter 24mm x 16mm



**90-1016- Marker Fixator**  
3.5mm



**90-1017 - Braunstien Marker**  
3.5 - 4.0mm




**90-3007 - Curette Double Ended**  
2mm & 3mm cup O/L130mm



**90-3008 - Curette Double Ended**  
2mm & 5mm cup O/L130mm



**90-3009 - Curette Double Ended**  
2mm & 4mm cup O/L130mm



90-8000 - Corneal Trephine 10mm



90-8005 - Corneal Trephine 15mm



90-8002 - Corneal Trephine 12mm



90-8006 - Corneal Trephine 16mm



90-8004 - Corneal Trephine 14mm



90-8007 - Corneal Trephine 17mm



91-6000 - Degree Gauge 11.7mm I.D.



91-6001 - Alignment Marker, 3, 6, 9 O'clock.



90-1020 - Thornton Ring



91-6002 - Toric IOL Marker 11mm



001 Sub Tenon Set




- 1 x 21-4009 Scissor Westcott Curved Rounded Tip
- 1 x 31-1006 Speculum Kratz-Barraquer Angled Temporal
- 1 x 10-1301 Forcep Suture Moorfields
- 1 x Wrap 50 x 50cm

001A Sub Tenon Set



- 1 x 21-4010 Scissor Westcott Extra Fine Curved Rounded Tip
- 1 x 31-1006 Speculum Kratz-Barraquer Angled Temporal
- 1 x 10-1301 Forcep Suture Moorfields
- 1 x Wrap 50 x 50cm

023 Intravitreoretinal Injection Set



- 1 x 30-1007 Speculum Kratz-Barraquer Straight
- 1 x 90-1015 Marker, Double Ended Scleral
- 1 x 10-1104 Forcep Straight Notched
- 1 x Wrap 50 x 50cm

208 Complete Ophthalmic Intravitreoretinal Set



- 1 x Tray with Integral Gallipot
- 1 x Trolley Cover
- 5 x Gauze 7.5 x 7.5cm
- 3 x Double Ended Cotton Buds
- 1 x Gallipot
- 1 x 90-1015 Double Ended Marker
- 1 x 31-2003 Speculum Barraquer Angled
- 1 x Plastic Rampley Sponger Holder
- 1 x 2684 Soft Touch Drape, Slit Incise, Opaque
- 3 x 9001 Eye Spears



092 - Chalazion Set



- 1 x 90-4005 Chalazion Clamp Small
- 1 x 90-3009 Curette Double Ended Med/Small
- 1 x 2333 No 11 Blade & Handle
- 1 x Wrap 50 x 50cm

111 - Chalazion Set



- 1 x 90-4006 Chalazion Clamp Medium
- 1 x 90-3009 Curette Double Ended Med/Small
- 1 x 2333 No.11 Blade on Handle
- 1 x Wrap 50 x 50cm

066 Suture Removal Set



- 1 x HR200 Scissor Iris
- 1 x HR125 Forcep Adson Non Toothed
- 1 x Wrap 50 x 55cm

038 Minor Dermatology Set



- 1 x Tray with Integral Gallipot
- 5 x Gauze 10 x 10cm
- 1 x Gallipot
- 1 x Hand Towel
- 1 x Trolley Cover
- 1 x HR300 BP Handle No.3
- 1 x HR519 Needle Holder Halsey Smooth
- 1 x 10-1302 Forcep Utility Lister 1:2 Teeth
- 1 x 20-7002 Scissor Iris Straight Sharp Tip
- 1 x 10-1006 Forcep Mosquito Straight
- 1 x 2467 Drape Aperture

151 Standard Cataract Set



- 1 x 41-5001 Manipulator Lens Hook Sinskey Angled
- 1 x 20-7002 Scissor Iris Straight Sharp Tip
- 1 x 10-1104 Forcep Suture Bonn 1:2 Teeth
- 1 x 31-1006 Speculum Kratz-Barraquer Angled Temporal
- 1 x 11-5024 Forcep Capsulorhexis Utrata
- 1 x Wrap 50 x 50cm

210 Ophthalmic Suture Removal Set



- 1 x 10-9007 Forcep Tying McPherson Straight
- 1 x 21-6004 Scissor Vannas Curved Blunt Tip
- 1 x Wrap 50 x 50cm

039 - Full Dermatology Set



- 1 x Tray with Integral Gallipot
- 5 x Gauze 10 x 10cm
- 1 x Gallipot
- 1 x Paper Towel
- 1 x Trolley Cover
- 1 x HR300 BP Handle No.3
- 1 x HR519 Needle Holder Halsey Smooth
- 1 x HR240 Scissor Curved Strabismus
- 1 x HR121 Forcep Adson Toothed
- 1 x HR125 Forcep Adson Non-Toothed
- 1 x HR250 Scissor Dressing Stitch Sharp/Blunt Tip
- 1 x HR940 Skin Hook Fine
- 1 x 20-7002 Scissor Iris Straight Sharp Tip
- 1 x 10-1006 Forcep Mosquito
- 1 x 10-1302 Forcep Utility Lister 1:2 Teeth
- 1 x 2467 Drape Aperture 50 x 55cm

047 - Minor Ops Set



- 1 x Tray with Integral Gallipot
- 10 x Gauze 10 x 10cm
- 2 x Gallipot
- 1 x 10-1302 Forcep Utility Lister
- 1 x 20-7002 Scissor Iris Straight Sharp Tip
- 1 x HR521 Needle Holder Halsey Smooth
- 1 x HR301 Handle Barron BP No.3
- 1 x HR400 Mosquito Curved
- 1 x HR122 Forcep Adson Toothed
- 1 x HR940 Skin Hook Fine
- 1 x HR283 Scissor Steven Tenotomy Curved Sharpe Tip
- 1 x 2467 Drape Aperture
- 1 x Wrap 50 x 50cm



074 - Dermatology Set



- 1 x Tray with Integral Gallipot
- 10 x Gauze 10 x 10cm
- 3 x Hand Towel
- 1 x Trolley Cover
- 1 x HR940 Skin Hook Fine
- 1 x HR301 Handle Barron BP No.3
- 1 x HR521 Needle Holder Halsey Smooth
- 1 x HR121 Forcep Adson Toothed
- 1 x HR283 Scissor Steven Tenotomy Curved Tip
- 1 x HR400 Mosquito Curved
- 1 x 2467 Drape Aperture

020 - Minor Ops Set



- 1 x Tray with Integral Gallipot
- 5 x Gauze 10 x 10cm
- 1 x Paper Towel
- 1 x HR281 Scissor Steven Tenotomy Straight
- 1 x HR200 Scissor Iris Straight Sharp Tip
- 1 x HR125 Forcep Adson
- 1 x HR121 Forcep Adson Toothed
- 1 x HR523 Needle Holder Nieverts
- 1 x HR301 Handle Barron BP No.3
- 1 x Plain Drape
- 1 x Wrap 50 x 55cm

- HR100 - Forcep Angled at 45° Jewellers**  
3mm tip O/L 110mm
- HR101 - Forcep Angled at 45° Jewellers**  
7mm tip O/L 110mm



- HR110 - Forcep Gillies Toothed**  
21mm Tying Platform O/L 110mm




- HR111 - Forcep De Bakey Toothed**  
20mm Tying Platform O/L 150mm



- HR120 - Forcep Adson Toothed**  
Dissecting 1.7mm Teeth O/L 120mm
- HR121 - Forcep Adson Toothed**  
Dissecting 0.8mm Teeth O/L 120mm
- HR122 - Forcep Adson Toothed**  
Dissecting 1.2mm Teeth O/L 120mm



- HR125 - Forcep Adson Non-Toothed**  
O/L 120mm



- HR126 - Forcep Adson Brown**  
O/L 120mm



- HR127 - Forcep Assistant Delicate**  
11mm Tying Platform O/L 115mm





**HR130 - Forcep Iris 11cm Serrated**  
15mm Tying Platform O/L 113mm



**HR131 - Forcep Epilation Angled**  
21mm Angled Tips O/L 91mm



**HR200 - Scissor Iris Straight**  
Sharp Tip 100mm O/L



**HR210 - Scissor Iris Curved**  
Sharp Tip O/L 100mm



**HR160 - Forcep Allis Tissue Tying 15cm**  
5:4 Teeth O/L 150mm



**HR150 - Forcep McIndoe 15cm**  
18mm Tying Platform O/L 155mm



**HR240 - Scissor Strabismus Curved**  
Rounded Tip O/L 115mm



**HR241 - Scissor Strabismus Straight**  
Rounded Tip O/L 115mm



**HR143 - Forcep Watch Maker No.7**  
13mm Tip O/L 115mm



**HR161 - Forcep Dissecting TOE 12.5cm**  
18mm Serrated Tip O/L 127zmm



**HR250 - Scissor Dressing Stitch**  
Sharp/Blunt Tip O/L 150mm



**HR261 - Scissor Mayo Curved**  
Rounded Tip O/L 185mm



**HR262 - Scissor Mayo Straight**  
Rounded Tip O/L 140mm



**HR270 - Scissor Metzemaum Curved 18cm (McIndoe)**  
Rounded Tip O/L 175mm



**HR300 - Handle Scalpel BP No.3**  
Rounded End O/L 125mm



**HR301 - Handle Barron BP No.3**  
O/L 135mm



**HR271 - Scissor Metzemaum Straight**  
Rounded Tip O/L 180mm



**HR280 - Scissor Straight**  
Sharp Tip O/L 125mm



**HR400 - Mosquito Curved**  
O/L 130mm  
**HR405 - Mosquito Curved**  
O/L 90mm



**HR410 - Mosquito Straight**  
O/L 90mm  
**HR420 - Mosquito Straight**  
O/L 130mm



## Mosquitos

**HR281 - Scissor Steven Tenotomy Straight**  
Sharp Tip O/L 110mm




**HR282 - Scissor Steven Tenotomy Curved**  
Supercut O/L 125mm



**HR283 - Scissor Steven Tenotomy Curved**  
Sharp Tip O/L 106mm



Should you have an idea for new single use surgical device, our in house design team would be happy to speak to you to see if we can turn it into reality. Please contact us at [sales@blinkmedical.com](mailto:sales@blinkmedical.com)





## Needle Holders

**HR500 - Needle Holder Mayo**  
O/L 150mm



**HR510 - Needle Holder Crilewood**  
O/L 130mm



The **Blink Medical** BiPolar range has now been increased to reflect the popularity of this instrument. The forceps which come complete with their leads offer a cost effective disposable alternative to the traditional reusable instrument. All forceps in the range are produced from the highest quality components and each one goes through multiple inspections here in the UK along with a comprehensive electrical test which ensures complete safety and efficacy.



**HR519 - Needle Holder Halsey**  
Smooth / Fine O/L 130mm



**HR520 - Needle Holder Halsey**  
O/L 130mm



**HR600 - Bipolar Mc Phearson**  
Straight with Lead O/L 110mm



**HR601 - Bipolar Mc Phearson**  
Straight with Lead O/L 200mm



**HR521 - Needle Holder Halsey**  
Smooth O/L 130mm



**HR522 - Needle Holder Kilner**  
O/L 130mm



**HR523 - Needle Holder Neiverts**  
O/L 130mm



**HR602 - Bipolar Adson Straight**  
with Lead O/L 150mm



**HR603 - Bipolar Bayonet**  
with Lead O/L 190mm



**HR700 - Retractor Catspaw Kilner**  
8mm Claw Diameter O/L 150mm



**HR703 - Dissector McDonald**  
O/L 195mm



**HR820 - Towel Clip Ball & Socket**  
O/L 100mm




**HR800 - Towel Clip Backhaus**  
O/L 95mm



**HR704 - Retractor self retaining**  
West (Sharp) O/L 135mm

**HR708 - Retractor self retaining**  
West (Blunt) 2:3 Teeth O/L 130mm



**HR776 - Hook Cawthorne**  
O/L 150mm



**HR951 - Sponge Holder Rampley**  
O/L 180mm



**HR777 - Retractor Catspaw Mini**  
5mm Claw Diameter O/L 165mm



**HR939 - St Bart's Hook**  
O/L 170mm



**HR940 - Skin Hook Fine**  
2mm Hook O/L 165mm

**HR941 - Skin Hook Standard**  
3mm Hook O/L 165mm



**HR830 - Curette Volkmann 5mm / 6mm cups**  
O/L 220mm





**HR608 - Cannula Fat Removal** V Shaped 9cm 1.5mm



**HR609 - Cannula Fat Injection** 9cm x 2mm



**HR610 - Cannula Fat Injection** 9cm x 1.5mm



**HR611 - Cannula Fat Injection** 12cm x 1.5mm



**HR612 - Cannula Fat Harvest** 15cm x 3mm



**HR613 - Cannula Infiltration** 30cm x 2.5mm, 10 ports, spiral pattern



**HR614 - Cannula Infiltration** 15cm x 2.5mm, 10 ports, spiral pattern



**HR615 - Cannula Infiltration** 20cm x 2mm, 10 ports, spiral pattern



**HR616 - Cannula Infiltration** 15cm x 2mm, 10 ports, spiral pattern



**HR617 - Cannula Infiltration** 10cm x 2mm, 10 ports, spiral pattern



All Blink Medical Lipo Cannulae are fitted with a Mercedes style tip and come in a variety of different lengths and diameters.

**HR620 - Cannula Liposuction** 30cm x 6mm



**HR621 - Cannula Liposuction** 30cm x 5mm



**HR622 - Cannula Liposuction** 30cm x 4mm



**HR623 - Cannula Liposuction** 30cm x 3mm



**HR624 - Cannula Liposuction** 15cm x 4mm



**HR625 - Cannula Liposuction** 15cm x 3mm



**HR626 - Cannula Liposuction** 15cm x 2mm



**HR643 - Cannula Liposuction** 22cm x 3mm - Luer Lock



**HR100 - Forcep Angled Jewellers**

3mm Tip O/L 110mm

**HR101 - Forcep Angled at 45° Jewellers**

7mm Tip O/L 110mm



**HR140 - Forcep Straight**

Serrated Tip 8mm O/L 95mm



**HR710 - Fine Croc Forcep**

6mm x 1.3mm Jaw



**HR711 - Forcep Tilley Nasal**

20mm Serrated Jaw



**HR141 - Forcep Angled Serrated Cutting**

Serrated Tip 7mm O/L 95mm



**HR142 - Forcep Jewellers No.5**

O/L 110mm



**HR714 - Scissor Crocodile**



**HR716 - Extra Fine Croc Forcep, Black Coating**

6mm x 1.2 jaw



**HR143 - Forcep Watch Maker No.7**

13mm Tip O/L 115mm



**HR715 - Extra Fine Croc Forcep**

5mm x 1.2mm Jaw



**HR720 - ENT Mirror**

25mm Mirror Diameter





**HR722 - Ear Speculum 2mm**

Plastic non sterile bag 100

**HR723 - Ear Speculum 3mm**

Plastic non sterile bag 100

**HR724 - Ear Speculum 4mm**

Plastic non sterile bag 100

**HR725 - Ear Speculum 5mm**

Plastic non sterile bag 100

**HR726 - Ear Speculum 6mm**

Plastic non sterile bag 100

**HR727 - Ear Speculum 7mm**

Plastic non sterile bag 100



**HR730 - Thudichum Small**

**HR731 - Thudichum Medium**

**HR732 - Thudichum Large**



**HR753 - Lempert Suction**

2mm O/D



**HR755 - Yankaure Suction Tube**



**HR740 - Tongue Depressor**



**HR750 - Zoellner Suction Tube**

Includes 20g Tip and Trocar



**HR756 - Fraiser Suction Tube**

4mm O/D



**HR761 - Forcep Troleitsch**

Small 5'



**HR751 - Jobson Horne Probe**

125mm O/L

**HR752 - Jobson Horne Probe 7"**

185mm O/L



**HR763 - Formby Cerumen Scoop and Hook**



**T013 - Ear Speculum 2mm**

**T014 - Ear Speculum 3mm**

**T015 - Ear Speculum 4mm**

**T016 - Ear Speculum 5mm**

**T017 - Ear Speculum 6mm**

**T018 - Ear Speculum 7mm**



**HR128 Baby Tischler Forcep**  
Shaft Length 210mm Jaw size 7 x 3 x 3mm



**HR133 - Forcep Eppendorfer Biopsy**  
Shaft length 200mm Jaw Size 7 x 3.5 x 4mm



**HR800 - Towel Clip Backhaus**



**HR181 - Forcep Bonney Polypus**



**HR134 - Forcep Tischler Morgan Biopsy**  
Shaft Length 205mm Jaw Size 6 x 5 x 3mm



**HR147 - Ultimate Eppendorfer**  
Shaft length 235mm Jaw Size 7 x 5 x 2mm



**P001 Nail Cutter**  
O/L 140mm Screw Joint



**P002 Nail Splitter**  
14cm Screw Joint



**HR830 - Curette Volkmann 5mm / 6mm cups**




**HR949 - Sponge Holder Rampley**  
240mm







2708 - 127 x 140cm Drape Blue Material with Fluid Collection pouch



2709 - 127 x 140cm Drape Opaque Material with Fluid Collection pouch



2707 - 52 x 26cm Drape Bi-Lateral Lasik, Daya with Fluid Collection pouch



Blink Medical are proud to announce the launch of their new ophthalmic drape range. Our Cataract Drapes feature a 360 degree fluid collection pouch rim, which ensures 100% capture of any operating fluid. This drape is available in Wetlaid or opaque plastic material.

The Daya Lasik Drape has been specially designed to facilitate improved and rapid bilateral patient draping during refractive and other corneal procedures. The drape has perforations in place avoiding the need for scissors and can replace an incise type drape. Additionally it features a small temporal pouch that collects fluid avoiding contamination and cleaning of the headrest and other areas between cases. Easy to use, pre-perforated and slightly larger than existing drapes, this is the idea drape for refractive surgery.

Page 1 - Introduction

Page 2 - Introduction

Page 3 - Company Information

Page 4 - Company Information

Page 5 - Ophthalmology Tissue Forceps

- 10-1415 - Forcep Birks Suture Tying
- 10-1302 - Forcep Utility Lister
- 10-1104 - Forcep Bonn
- 10-1105 - Forcep Bonn
- 10-9007 - Forcep Tying McPherson
- 11-7009 - Forcep Piers Colibri Notched
- 11-7010 - Forcep Colibri Toothed
- 11-7011 - Forcep Colibri Toothed

Page 6 - Ophthalmology Tissue Forceps

- 10-1203 - Forcep St Martins
- 10-1204 - Forcep St Martins
- 10-1205 - Forcep Barraquer Fixation
- 10-9008 - Forcep Tying McPherson
- 10-8008 - Forcep Hoskins Straight Notched
- 11-6010 - Forcep Tying Kelman McPherson

Page 7 - Ophthalmology Capsulorhexis Forceps

- 11-5011 - Forcep Capsulorhexis Inamura
- 11-5012 - Forcep Capsulorhexis Utrata
- 11-5021 - Forcep Capsulorhexis Inamura
- 11-5022 - Forcep Capsulorhexis Utrata
- 11-5025 - Forcep Capsulorhexis Fulcher

Page 8 - Ophthalmology Capsulorhexis Forceps

- 11-5013 - Forcep Capsulorhexis Utrata
- 11-5020 - Forcep Capsulorhexis Inamura
- 11-5023 - Forcep Capsulorhexis Micro Incision 1.8mm
- 11-5024 - Forcep Capsulorhexis Utrata Curved

Page 9 - Ophthalmology Forceps Other

- 10-1300 - Forcep Epilation
- 10-1301 - Forcep Suture Moorfields
- 10-1311 - Forcep Suture Moorfields Utility
- 11-1005 - Forcep Mosquito Curved
- 11-1006 - Forcep Mosquito Straight

Page 10 - Ophthalmology Scissors

- 20-7001 - Scissor Iris Straight
- 20-7002 - Scissor Iris Straight
- 20-6003 - Scissor Vannas Straight
- 21-6004 - Scissor Vannas Curved
- 20-4006 - Scissor Westcott Straight Blunt
- 20-4007 - Scissor Westcott Straight Sharp

Page 11 - Ophthalmology Scissors

- 21-4008 - Scissor Westcott Curved Sharp
- 21-4009 - Scissor Westcott Curved Blunt
- 21-4010 - Scissor Westcott Curved Blunt
- 21-5001 - Scissor Westcott Curved Small
- 21-3010 - Scissor Ong

Page 12 - Ophthalmology Lens Chopping & Chopping

- 40-5002 - Manipulator Lens Hook Sinskey Straight
- 41-5001 - Manipulator Lens Hook Sinskey Angled
- 41-3004 - Manipulator Spatula
- 41-2005 - Rotator Nucleus Y Shaped
- 41-1006 - Lester IOL Manipulator
- 41-7008 - Chopper Greens Phaco
- 41-7009 - Chopper Claw
- 41-4005 - Chopper Sadiq
- 41-4003 - Manipulator Phaco Chopper Mackool
- 41-3005 - Neuhaus Spatula
- 41-2006 - Rotator Nucleus Y Shaped
- 41-6007 - Paddle Manipulator Drysdale
- 41-8009 - Manipulator Mushroom Tipped
- 41-4004 - Chopper Ambidextrous
- 41-4006 - Spatula Govan
- 41-4007 - Haeffiger Chopper Sadiq
- 41-9010 - Manipulator Melonies Double Ended

Page 13 - Ophthalmology Speculae

- 30-4000 - Speculum Barraquer
- 30-4001 - Speculum Adjustable
- 30-3001 - Speculum Burratto
- 30-3002 - Speculum Paediatric Spring
- 31-3003 - Speculum Cook Adjustable
- 31-3004 - Speculum Clarkes
- 31-3100 - Speculum Neo Natal
- 31-3101 - Speculum Neo Natal
- 31-3002 - Speculum Leibinger

Page 14 - Ophthalmology Speculae

- 31-2003 - Speculum Barraquer
- 30-2004 - Speculum Barraquer Straight
- 31-1006 - Speculum Kratz-Barraquer
- 30-1007 - Speculum Kratz-Barraquer

Page 15 - Ophthalmology Refractive / Vitreo

- 51-0000 - Flap Lifter Lasik
- 51-0001 - Epithelial Spatula
- 51-0002 - Alcohol Well

Page 16 - Ophthalmology Vitreoretinal

- 60-1000 - Forcep End Gripping
- 60-1001 - Forcep ILM
- 60-1002 - A-Symmetrical
- 60-1003 - Forcep Straight
- 60-2000 - Scissor Straight

Page 17 - Ophthalmology Oculoplastic

- 90-2008 - Probe Liebricht 000+00
- 90-2009 - Probe Liebricht 00+0
- 90-2010 - Probe Liebricht 1+2
- 90-2011 - Probe Liebricht 3+4
- 90-1011 - Dilator Punctal Double Ended
- 90-1012 - Nettleships Dilator
- 90-8012 - Bulldog Clamp

Page 18 - Ophthalmology Oculoplastic

- 90-7001 - Needle Holder Castroviejo
- 90-7002 - Needle Holder Castroviejo
- 91-7003 - Needle Holder Castroviejo
- 90-7004 - Needle Holder Troutman
- 90-5004 - Hook Stabismus Graefe Medium
- 90-5005 - Hook Alexandra Muscle

Page 19 - Ophthalmology Oculoplastic

- 90-4005 - Clamp Chalazion Small
- 90-4006 - Clamp Chalazion Medium
- 90-4007 - Clamp Chalazion Large
- 90-4008 - Desmares Entropion Clamp Left
- 90-4009 - Desmares Entropion Clamp Right
- 90-3007 - Curette Double Ended
- 90-3008 - Curette Double Ended
- 90-3009 - Curette Double Ended

Page 20 - Ophthalmology for Intravitreal Injections

- 90-1014 - Caliper and Rule
- 90-1015 - Double Ended Marker
- 90-1016 - Marker Fixator
- 90-1017 - Braunstien Marker

Page 21 - Ophthalmology Other

- 90-9000 - Corneal Trephine 10mm
- 90-8002 - Corneal Trephine 12mm
- 90-8004 - Corneal Trephine 14mm
- 90-8005 - Corneal Trephine 15mm
- 90-8006 - Corneal Trephine 16mm
- 90-8007 - Corneal Trephine 17m
- 91-6000 - Degree Gauge 11.7mm I.D.
- 91-6001 - Alignment Marker, 3, 6, 9 O'clock.
- 90-1020 - Thornton Ring
- 91-6002 - Toric IOL Marker 11mm

Page 22 - Ophthalmology Standard Ophthalmic Sets

- 001 - Sub Tenon Set
- 001A - Sub Tenon Set
- 023 - Intravitreal Injection Set
- 208 - Complete Ophthalmic Intravitreal Set

Page 23 - Ophthalmology Standard Ophthalmic Sets

- 092 - Chalazion Set
- 111 - Chalazion Set
- 151 - Standard Cataract Set
- 210 - Ophthalmic Suture Removal Set

Page 24 - General Surgery Standard Sets

- 066 - Suture Removal Set
- 036 - Minor Dermatology Set
- 039 - Full Dermatology Set
- 047 - Minor Ops Set

Page 25 - General Surgery Standard Sets

- 074 - Dermatology Set
- 020 - Minor Ops Set

Page 26 - General Surgery Forceps

- HR100 - Forcep Angled at 45° Jewellers
- HR723 - Ear Speculum 3mm
- HR101 - Forcep Angled at 45° Jewellers
- HR110 - Forcep Gillies Toothed
- HR111 - Forcep De Bakey Toothed
- HR120 - Forcep Adson Toothed
- HR121 - Forcep Adson Toothed
- HR122 - Forcep Adson Toothed
- HR125 - Forcep Adson Non-Toothed
- HR126 - Forcep Adson Brown
- HR127 - Forcep Assistant Delicate

Page 27 - General Surgery Forceps

- HR130 - Forcep Iris 11cm Serrated
- HR131 - Forcep Epilation Angled
- HR160 - Forcep Allis Tissue Tying 15cm
- HR150 - Forcep McIndoe 15cm
- HR143 - Forcep Watch Maker No.7
- HR161 - Forcep Dissecting TOE 12.5cm

Page 28 - General Surgery Scissors

- HR200 - Scissor Iris Sharp Straight
- HR210 - Scissor Iris Sharp Curved
- HR240 - Scissor Curved Strabismus
- HR241 - Scissor Strabismus Straight
- HR250 - Scissor Dressing Stitch Sharp/Blunt
- HR261 - Scissor Mayo Curved

Page 29 - General Surgery Scissors

- HR262 - Scissor Mayo Straight
- HR270 - Scissor Metzemabaum Curved 18cm (McIndoe)
- HR271 - Scissor Metzemabaum Straight
- HR280 - Scissor Sharp Straight
- HR281 - Scissor Steven Tenotomy
- HR282 - Scissor Steven Tenotomy Curved
- HR283 - Scissor Steven Tenotomy Curved

Page 30 - General Surgery Handles / Mosquitos

- HR300 - Handle Scalpel BP No.3
- HR301 - Handle Barren BP No.3
- HR400 - Mosquito Curved
- HR405 - Mosquito Curved
- HR410 - Mosquito Straight
- HR420 - Mosquito Straight

Page 31 - General Surgery Needle Holders

- HR500 - Needle Holder Mayo
- HR510 - Needle Holder Crilewood
- HR519 - Needle Holder Halsey
- HR520 - Needle Holder Halsey
- HR521 - Needle Holder Halsey
- HR522 - Needle Holder Kilner
- HR523 - Needle Holder Neiverts

Page 32 - General Surgery Bipolar

- HR600 - Bipolar Mc Phearson
- HR601 - Bipolar Mc Phearson
- HR602 - Bipolar Adson Straight
- HR603 - Bipolar Bayonet

Page 33 - General Surgery Retractors

- HR700 - Retractor Catspaw Kilner
- HR703 - Dissector McDonald
- HR704 - Retractor self retaining
- HR708 - Retractor self retaining
- HR776 - Hook Cawthorne
- HR777 - Retractor Catspaw Mini
- HR939 - St Bart's Hook
- HR940 - Skin Hook Fine
- HR941 - Skin Hook Standard

Page 34 - - General Surgery Retractors

- HR820 - Towel Clip Ball & Socket
- HR800 - Towel Clip Backhaus
- HR951 - Sponge Holder Rampley
- HR830 - Curette Volkmann 5mm / 6mm cups

Page 35 - Cosmetic Surgery Cannulae

- HR608 - Cannula Fat Removal
- HR609 - Cannula Fat Injection
- HR610 - Cannula Fat Injection
- HR611 - Cannula Fat Injection
- HR612 - Cannula Fat Harvest
- HR613 - Cannula Infiltration
- HR614 - Cannula Infiltration
- HR615 - Cannula Infiltration
- HR616 - Cannula Infiltration
- HR617 - Cannula Infiltration

Page 36 - Cosmetic Surgery Cannulae

- HR620 - Cannula Liposuction
- HR621 - Cannula Liposuction
- HR622 - Cannula Liposuction
- HR623 - Cannula Liposuction
- HR624 - Cannula Liposuction
- HR625 - Cannula Liposuction
- HR626 - Cannula Liposuction
- HR643 - Cannula Liposuction

Page 37 - Cosmetic Surgery Forceps

- HR100 - Forcep Angled Jewellers
- HR101 - Forcep Angled at 45° Jewellers
- HR140 - Forcep Straight
- HR141 - Forcep Angled Serrated Cutting
- HR142 - Forcep Jewellers No.5
- HR143 - Forcep Watch Maker No.7

Page 38 - ENT

- HR710 - Fine Croc Forcep
- HR711 - Forcep Tilley Nasal
- HR714 - Scissor Crocodile
- HR716 - Extra Fine Croc Forcep, Black Coating
- HR715 - Extra Fine Croc Forcep
- HR720 - ENT Mirror

Page 39 - ENT

- HR722 - Ear Speculum 2mm
- HR723 - Ear Speculum 3mm
- HR724 - Ear Speculum 4mm
- HR725 - Ear Speculum 5mm
- HR726 - Ear Speculum 6mm
- HR727 - Ear Speculum 7mm
- HR730 - Thudichum Small
- HR731 - Thudichum Medium
- HR732 - Thudichum Large
- HR740 - Tongue Depressor

Page 40 - ENT

- HR753 - Lempert Suction Retracted
- HR755 - Yankauer Suction Tube Retracted
- HR756 - Fraiser Suction Tube
- HR761 - Forcep Trolletsch
- HR763 - Cerumen Scoop
- T013 - Ear Speculum 2mm
- T014 - Ear Speculum 3mm
- T015 - Ear Speculum 4mm
- T016 - Ear Speculum 5mm
- T017 - Ear Speculum 6mm
- T018 - Ear Speculum 7mm

Page 41 - Other Other

- HR128 - Baby Tichler Forcep
- HR133 - Forcep Eppendorfer Biopsy
- HR134 - Forcep Tichler Morgan Biopsy
- HR147 - Ultimate Eppendorfer
- HR830 - Curette Volkmann 5mm / 6mm cups

Page 42 - Other Other

- HR800 - Towel Clip Backhaus
- HR181 - Forcep Bonney Polypus
- P001 - Nail Cutter
- P002 - Nail Splitter
- HR949 - Sponge Holder Rampley

Page 43 - Drapes

- 2708 - 127 x 140cm Drape Blue Material with Fluid Collection pouch
- 2709 - 127 x 140cm Drape Opaque Material with Fluid Collection pouch
- 2707 - 52 x 26cm Drape Bi-Lateral Lasik, Daya with Fluid Collection pouch

Page 44 - Index

This image shows a single page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page, leaving small margins at the top and bottom. There are no vertical margin lines, and the page is completely blank except for the lines themselves.